

Re-Packaging your ForceFrame Fold

Once you have received your ForceFrame Max, you will need to re-purpose the packaging for sending back your ForceFrame Fold.


You will need:

- ForceFrame Max box (Box A) & Crossbar box (Box B) & inner packing
- 5mm Allen key
- 13mm spanner
- Stanley knife
- Packing tape

Disassembling your ForceFrame Fold – Box A

Note: Ensure your Head Unit is turned off for transit


1. Remove the Top Bar by unscrewing the wing bolts by hand – screw the bolts back together and place them with the Top Bar to the side




2. Remove the Crossbar by pulling both the vertical indexing plungers, either side of the Crossbar, outward and lifting the Crossbar up and off the top of the Posts. Place the Crossbar to the side




3. Remove 2 x Frame Locking Plates by unscrewing the 8 x wing bolts, by hand from either side of the ForceFrame. Place Frame Locking Plates and wing bolts to the side



4. Remove the 2 x Posts from the Base (4 bolts per Post) using 5mm Allen key . Attach the bolts back onto the centre Base piece once they are removed. Place both Posts to the side



5. Separate the 3 pieces of the ForceFrame Base (2 bolts per side) using 5mm allen key and 13mm spanner. Attach bolts back into the centre base piece



Packing your ForceFrame Fold into your ForceFrame Max Box



1. Reuse 2 x foam inserts to pad the bottom of the ForceFrame Max box (Box A). Using a Stanley knife, cut them through the centre as per the photos below.

Note: It is essential to cut these, as not cutting these will not allow enough room to close the box at the end of repacking.



2. Place the 4 x cut pieces of foam into the base of the box (Box A)



3. Place the large Base Plate into the box (Box A) against one side



4. Place the smaller Base Plate into the box (Box A) against the large base plate




5. Place the additional 2 pieces of foam to each side of the box (Box A).




6. Place the centre join of the Base Plate into the box (Box A) and push down into the slits of the foam (to hold it in place).



7. Tape the 2 wing bolts to the top bar, and place this into the box (Box A) on top of the centre piece of the Base Plate.



8. Place the 8 x wing bolts onto small piece of bubble wrap and fold up. Place 2 x Frame Locking Plates onto bubble wrap, place the wing bolts on top and fold both together. Place into the box (Box A)





9. Tape the 2 x Posts together and place on top of the Top Bar in the box (Box A)



10. Place 2 x corner foam pieces on either side of the top corners of the Base Plate. Push down gently to ensure they are flush with the top of the box (Box A). Place any additional bubble wrap into the additional space to protect the base plate and other parts.




11. Tape the box up (Box A) securely with tape



Packing your Crossbar – Box B

1. Use 2 x foam inserts to place on each end of the Crossbar



2. Insert Crossbar into the box (Box B) and tape the box up securely. Tape the box up (Box B) securely with tape



3. Contact support@vald.com to initiate how to send your Crossbar and ForceFrame Fold back to VALD.